


Nitinol Flexible Reamers

MODULAR SYSTEM

Modular Drive Shaft \varnothing 5 mm - TL 470 mm

Head sizes \varnothing 9 mm - \varnothing 19.5 mm in 1/2 mm increments
(Reamer heads to \varnothing 24 mm available on request)

Guide Wire \varnothing 3 mm, CoCr or SS



Power Tool Fittings

AO
Stryker/Zimmer-Hall
Hudson

MONOBLOC SYSTEM

for femoral application

Monobloc Reamers \varnothing 4.95 mm
Total Length 475 mm
Head sizes \varnothing 8-13 mm in 1/2 mm increments
Guide Wire \varnothing 3 mm, CoCr or SS

for humeral application

Monobloc Reamers \varnothing 4 mm
Total Length 350 mm
Sizes \varnothing 6-10 mm in 1/2 mm increments
(\varnothing 10.5 mm and + available on request)
Guide Wire \varnothing 2 mm, CoCr or SS