



Designed to **facilitate the conversion from temporary to definitive fixation**

Designed for **a minimally invasive approach**

Compatible with Galaxy Fixation System

Designed to **avoid contamination of the tibial medullary canal**

Kickstand
for foot elevation option



Easy Closure
by Allen Wrench 5mm

Radiolucent



All sterile components can be ordered separately or grouped by procedure through the **Galaxy UNYCO Boxes**



99-93574
Galaxy UNYCO
Mini Kit Tibia Sterile

- 1x Large Multiscrew Clamp
- 4x UNYCO Cancellous Screw



99-93509
Galaxy UNYCO
Mini Kit Instruments Sterile

- 1x Allen Wrench
- 1x Power Drill Torque Limiter



99-93477
Galaxy UNYCO
Mini Kit Ankle Unilateral Sterile

- 1x Radiolucent Foot Unit
- 3x UNYCO Cancellous Screw



99-93570
Galaxy UNYCO
Cancellous Screw Kit Sterile

- 2x UNYCO Cancellous Screw

99-93794
Galaxy UNYCO
Diaphyseal Tibia Box

- 2x Mini Kit Tibia Sterile
- 1x Mini Kit Instruments Sterile
- 1x Galaxy Rod Ø 12mm L 350mm sterile

99-93793
Galaxy UNYCO
Ankle Bridging Box -
Delta Frame


- 1x Mini Kit Tibia Sterile
- 1x Mini Kit Instruments Sterile
- 2x Transitional Clamp Sterile
- 1x Large Single Clamp Sterile
- 2x Galaxy Rod Ø 12mm L 350mm sterile
- 1x Transfixing Pin Ø 4mm Sterile

99-93792
Galaxy UNYCO
Ankle Bridging Box -
Unilateral Frame

- 1x Mini Kit Tibia Sterile
- 1x Mini Kit Ankle Unilateral Sterile
- 1x Mini Kit Instruments Sterile
- 6x Large Single Clamp Sterile
- 2x Galaxy Rod Ø 12mm L 350mm sterile



Distributed by:



Manufactured by:
ORTHOFIX Srl
Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy
Telephone +39 045 6719000, Fax +39 045 6719380



@galaxyunyc
www.galaxyunyc.com

www.orthofix.com

MC-1511-PL-E0 A 10/15

Instructions for Use: See actual package insert for Instructions for Use.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician. Proper surgical procedure is the responsibility of the medical professional. Operative techniques are furnished as an informative guideline. Each surgeon must evaluate the appropriateness of a technique based on his or her personal medical credentials and experience. Please refer to the "Instructions for Use" supplied with the product for specific information on indications for use, contraindications, warnings, precautions, adverse reactions and sterilization.

