



MINIFIXATION



MultiPlanar MiniRail™



Multiplanar
Adjustment


Mini Fixation

Micrometric
Correction



MultiPlanar MiniRail™


The geared, M2 MultiPlanar MiniRail permits staged correction of complex joint deformities and non-surgical lengthening of musculotendinous contractures. With the M2 positioned over the center of rotation of the deformity, gradual correction in both the radio-ulnar and dorso-palmar planes and compression/distraction of either railed segment is possible.



- Multiplanar Adjustment
- Micrometric Correction
- Mini Fixation

- Geared rails allow gradual positioning in both radio-ulnar or dorso-palmar planes
- Distraction and compression can be accomplished in either railed segment

Radial Club Hand Application



CODE	INDICATIONS
M511	<ul style="list-style-type: none"> • Complex upper and lower extremity deformities • Hand, wrist or foot contractures • Indicated for pediatrics and adult patients

Adult & Pediatric Deformity Correction | Foot & Ankle | Bone Growth Stimulation | Fixation | Biologics