



M<sup>2</sup> MULTIPLANAR MINIRAIL



MultiPlanar MiniRail™



Multiplanar  
Adjustment

Mini Fixation

Micrometric  
Correction




ORTHOFIX®



# MultiPlanar MiniRail™


The geared, M2 MultiPlanar MiniRail permits staged correction of complex joint deformities and non-surgical lengthening of musculotendinous contractures. With the M2 positioned over the center of rotation of the deformity, gradual correction in both the radio-ulnar and dorso-palmar planes and compression/distraction of either railed segment are possible.



- **Multiplanar Adjustment**
- **Micrometric Correction**
- **Mini Fixation**

- Geared rails allow **gradual** positioning in both radio-ulnar or dorso-palmar planes
- Distraction and compression can be accomplished in either railed segment

Radial Club Hand Application



CODE	INDICATIONS
M511	<ul style="list-style-type: none"> <li>• Complex upper and lower extremity deformities</li> <li>• Hand, wrist or foot contractures</li> <li>• Indicated for paediatrics and adult patients</li> </ul>

Deformity Correction | Trauma | Pediatrics | Bone Growth Stimulation



Your Distributor is: