

RadioLucent Wrist Fixator

Convenience
at hand



Excellent fracture visualisation, unobstructed view of radius and carpals



Ball-joint design

helps anatomical reduction

Low profile design

for greater patient comfort

Compression-distraction unit

for accurate reduction



Telescoping radiolucent body

Ready to use fixator

in convenient sterile kit including required instrumentation, compression-distraction unit and standard bone screws



Your Distributor is:



Manufactured by: ORTHOFIX Srl.
Via Delle Nazioni 9 - 37012 Bussolengo (Verona) - Italy
Telephone 0039-0456719000 - Fax 0039-0456719380

Ordering Information

Code No	Description
99-36501	RadioLucent Wrist Fixator Kit, sterile
For further ordering information see section "RadioLucent Wrist Fixator" in the Orthofix General Catalogue	

See Quick Reference Guide PG 150

References:

1) Nelson D.L.: "The RadioLucent Wrist Fixator for Distal Radius Fractures" Orthofix External Fixation in Trauma and Orthopaedics, G. De Bastiani et al. Springer 2000, Ch. 17: 181-191.